



Recombinant Arabidopsis thaliana F-box/LRR-repeat protein 13 (FBL13)

| | |
|--------------------------|--|
| Product Code | CSB-MP880837DOA |
| Abbreviation | FBL13 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9FWZ1 |
| Product Type | Recombinant Protein |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MDEDREKHVR AKGSDEVDWI SKLPDCLLCE VLLNLPTKDV VKTSVLSRRW RNLWKHVPGL DLDNTDFQEF NTFLSFVDSF LDFNSESFLQ KFILKYDCDD EYDPDIFLIG RWINTIVTRK VQHIDVLDD S YGSWEVQLPS SIYTCELSVS LKLCGLTLAS PEFVSLPSLK VMDLIITKFA DDMGLETIT KCPVLESLTI ERSFCDEIEV LRVRSQSLR FTHVADSDEG VVEDLVV SID APKLEYLRLS DHRVASFILN KPGKLVKADI DIVFNLSVN KFNPDLPKR TMIRNLLGI STIKDMIIFS STLEVIYDFS RCERLPLFRN LSVLCVEFYG YMWEMLP IFL ESCPNLKTLV VKSASYQEK G ENIILPGPRR FLSSLEYVKI ERPLKGEAME MKLVSYLLEN STILKKLTL C LDDSVKKEDS VILKELLAIP RLSTSSKVVV L |
| Source | Mammalian cell |
| Target Names | FBL13 |
| Protein Names | Recommended name: F-box/LRR-repeat protein 13 |
| Expression Region | 1-451 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |