



Recombinant Arabidopsis thaliana Germin-like protein subfamily 2 member 2 (At1g02335)

Product Code	CSB-BP880852DOA
Abbreviation	At1g02335
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FZ27
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	YDPDALQD LCVADKSHGT KLNGFPCKET LNITESDFFF AGISKPAVIN STMGSAVTGA NVEKIPGLNT LSVSLARIDY APGGLNPPHT HPRATEVVYV LEGELEVGF I TTANKLFTKT IKIGEVFVFP RGLVHFQKNN GKSPASVLSA FNSQLPGTAS VAATLFAAEP ALPEDVLTKT FQVGSKMVDK IKERLATKK
Source	Baculovirus
Target Names	At1g02335
Protein Names	Recommended name: Germin-like protein subfamily 2 member 2
Expression Region	23-219
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.