



Recombinant Galago senegalensis Somatotropin (GH1)

Product Code	CSB-YP880885GAB
Abbreviation	GH1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9GKA1
Product Type	Recombinant Protein
Immunogen Species	Galago senegalensis (Northern lesser bushbaby) (Senegal bushbaby)
Purity	≥85% (SDS-PAGE)
Sequence	AFPA MPLSSLFANA VLRAQHLHQL AADTYKEFER AYIPEGQRYS IQNTQAAFCF SETIPAPTKG DEAQQRSDME LLRFSLLLIQ SWLGPVQLLS RVFTNSLVLG TLDRVYEKLG DLEEGIQALM RELEDGSPRV GQILKQTYDK FDTNLRSDDA LLKNYGLLSC FKKDLHKAET YLRVMKCRRF VESSCAF
Source	Yeast
Target Names	GH1
Protein Names	Recommended name: Somatotropin Alternative name(s): Growth hormone
Expression Region	27-217
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.