



Recombinant Human Brain-specific serine protease 4 (PRSS22)

Product Code	CSB-MP880923HU
Abbreviation	PRSS22
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9GZN4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ARIPVPPA CGKPQQLNRV VGGEDSTDSE WPWIVSIQKN GTHHCAGSLL TSRWVITAAH CFKDNLNKPY LFSVLLGAWQ LGNPGSRSQK VGVAVVEPH VYSWKEGACA DIALVRLERS IQFSERVLPI CLPDASIHLP PNTHCWISGW GSIQDGVPLP HPQTLQKLKV PIIDSEVCSH LYWRGAGQGP ITEDMLCAGY LEGERDACLG DSGGPLMCQV DGAWLLAGII SWGEGCAERN RPGVYISLSA HRSWVEKIVQ GVQLRGRAQG GGALRAPSQG SGAAARS
Source	Mammalian cell
Target Names	PRSS22
Protein Names	Recommended name: Brain-specific serine protease 4 Short name= BSSP-4 EC= 3.4.21.- Alternative name(s): Serine protease 22 Serine protease 26 Trypsin epsilon
Expression Region	33-317
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a member of the trypsin family of serine proteases. The enzyme is expressed in the airways in a developmentally regulated manner. The gene is part of a cluster of serine protease genes on chromosome 16.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.