



Recombinant Human NIF3-like protein 1 (NIF3L1)

Product Code	CSB-MP880926HU
Abbreviation	NIF3L1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9GZT8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLSSCVRPVP TTVRFVDSL I CNSSRSFMDL KALLSSLNDF ASLSFAESWD NVGLLVEPSP PHTVNTLFLT NDLTEEVME E VLQKKADLIL SYHPPIFRPM KRITWNTWKE RLVIRALENR VGIYSPHTAY DAAPQGVNNW LAKGLGACTS RPIHPSKAPN YPTEGNHRVE FNVNYTQDLD KVMSAVKGID GVSVTSFSAR TGNEEQTRIN LNCTQKALMQ VVDFLSRNKQ LYQKTEILSL EKPLLLHTGM GRLCTLDESV SLATMIDRIK RHLKLSHIRL ALGVGRTLES QVKVVALCAG SGSSVLQGVE ADLYLTGEMS HHDTLDAASQ GINVILCEHS NTERGFLSDL RDMLDSHLEN KINIILSETD RDPLQVV
Source	Mammalian cell
Target Names	NIF3L1
Protein Names	Recommended name: NIF3-like protein 1 Alternative name(s): Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 1 protein
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.