



Recombinant Human N-alpha-acetyltransferase 50 (NAA50)

Product Code	CSB-BP880928HU
Abbreviation	NAA50
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9GZZ1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MKGSRIELGD VTPHNIKQLK RLNQVIFPVS YNDKFYKDVLEVGELAKLAY FNDIAVGAVC CRVDHSQNQK RLYIMTLGCL APYRRLGIGT KMLNHVLNIC EKDGTFDNIY LHVQISNESA IDFYRKFGFE IETKKNYYK RIEPADAHVL QKNLKVPSGQ NADVQKTDN
Source	Baculovirus
Target Names	NAA50
Protein Names	Recommended name: N-alpha-acetyltransferase 50 EC= 2.3.1.- Alternative name(s): N-acetyltransferase 13 N-acetyltransferase 5 Short name= hNAT5 N-acetyltransferase san homolog Short name= hSAN NatE catalytic subunit
Expression Region	1-169
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.