



Recombinant Human Coiled-coil domain-containing protein 113 (CCDC113)

Product Code	CSB-YP880934HU
Abbreviation	CCDC113
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H0I3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTDDESESVL SDSHEGSELE LPVIQLCGLV EELSYVNSAL KTETEMFEKY YAKLEPRDQR PPRLSEIKIS AADYAQFRGR RRSKSRTGMD RGVGLTADQK LELVQKEVAD MKDDL RHTRA NAERDLQHHE AIIEEAEIRW SEVSREVEHF EKDILKAISK KKGSI L ATQK VMKYIEDMNR RRDNMKEKLR LKNVSLKVQR KKMLLQLRQK EEVSEALHDV DFQQLKIENA QFLETI EARN QELTQLKLSS GNTLQVLNAY KSKLHKAMEI YLNLDKEILL RKELLEKIEK ETLQVEEDRA KAEAVNKRLR KQLAEFRAPQ VMTYVREKIL NADLEKSIRM WERKVEIAEM SLKGHRKAWN RMKITNEQLQ ADYLAGK
Source	Yeast
Target Names	CCDC113
Protein Names	Recommended name: Coiled-coil domain-containing protein 113
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.