



Recombinant Human Coiled-coil domain-containing protein 8 (CCDC8)

Product Code	CSB-EP880938HU-B
Abbreviation	CCDC8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H0W5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLQIGEDVDY LLIPREVRLA GGVWRVISKP ATKEAEFRER LTQFLEEEGR TLEDVARIME KSTPHPPQPP KKPKEPRVRR RVQQMVTPPP RLVVGTYDSS NASDSEFSDF ETSRDKSRQG PRRGKKVRKM PVSYLGSKFL GSDLESEDDE ELVEAFLRRQ EKQPSAPPAR RRVNLPVPMF EDNLGPQLSK ADRWREYVSQ VSWGKLRKRRV KGWAPRAGPG VGEARLASTA VESAGVSSAP EGTSPGDRLG NAGDVCVPQA SPRRWRPKIN WASFRRRRKE QTAPTGGAD IEADQGGEEA DSQREEAIAD QREGAAGNQR AGAPADQGAE AADNQREEAA DNQRAGAPAE EGAEAADNQR EEAADNQRAE APADQRSQGT DNHREEAADN QRAEAPADQG SEVTDNQREE AVHDQRRERAP AVQGADNQRA QARAGQRAEA AHNQRAGAPG IQEAEVSAAQ GTTGTAPGAR ARKQVKTVRF QTPGRFSWFC KRRRAFVWHTP RLPTLPKRVP RAGEARNLRV LRAEARAEAE QGEQEDQL
Source	E.coli
Target Names	CCDC8
Protein Names	Recommended name: Coiled-coil domain-containing protein 8
Expression Region	1-538
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.