



Recombinant Human Melanoma-associated antigen H1 (MAGEH1)

Product Code	CSB-MP880949HU
Abbreviation	MAGEH1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H213
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPRGRKSRRR RNARAAEENR NNRKIQASEA SETPMAASVV ASTPEDDLTG PEEDPSTPEE ASTTPEEASS TAQAQKPSVP RSNFQGTTKS LLMSILALIF IMGNSAKEAL VWKVLGKLG M QPGRQHSIFG DPKKIVTEEF VRRGYLIYKP VPRSSPVEYE FFWGPRAHVE SSKLKVMHFV ARVRNRCSKD WPCNYDWDSD DDAEVEAILN SGARGYSAP
Source	Mammalian cell
Target Names	MAGEH1
Protein Names	Recommended name: Melanoma-associated antigen H1 Alternative name(s): Apoptosis-related protein 1 Short name= APR-1 MAGE-H1 antigen Restin
Expression Region	1-219
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene is thought to be involved in apoptosis. Multiple polyadenylation sites have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.