



Recombinant Human Coiled-coil domain-containing protein 68 (CCDC68)

Product Code	CSB-EP880954HU-B
Abbreviation	CCDC68
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H2F9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTTVTVTTEI PPRDKMEDNS ALYESTSAHI IEETEYVKKI RTTLQKIRTQ MFKDEIRHDS TNHKLDAKHC GNLQQGSDSE MDPSCCSLDL LMKKIKGKDL QLLEMNKENE VLKIKLQASR EAGAAALRNV AQLRFENYQT QSEEVRRKKQE DSKQLLQVNK LEKEQKLKQH VENLNQVAEK LEEKHSQITE LENLVQRMEK EKRTLLERKL SLENKLLQLK SSATYGKSCQ DLQREISILQ EQISHLQFVI HSQHQNLRSV IQEMEGLKNN LKEQDKRIEN LREKVNILEA QNKELKTQVA LSSETPRTKV SKAVSTSELK TEGVSPYLML IRLRK
Source	E.coli
Target Names	CCDC68
Protein Names	Recommended name: Coiled-coil domain-containing protein 68 Alternative name(s): Cutaneous T-cell lymphoma-associated antigen se57-1 Short name=CTCL-associated antigen se57-1
Expression Region	1-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.