



Recombinant Human Ubiquilin-3 (UBQLN3)

Product Code	CSB-MP880963HU
Abbreviation	UBQLN3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H347
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAKGGEALPQ GSPAPVQDPH LIKVTVKTPK DKEDFSVTDT CTIQQLKEEI SQRFKAHPDQ LVLIFAGKIL KDPDSLAQCG VRDGLTVHLV IKRQHRAMGN ECPAASVPTQ GSPGSLPQP SSIYPADGPP AFSLGLLTGL SRLGLAYRGF PDQPSSLMRQ HVSVPEFVTQ LIDDPFIPGL LSNTGLVRQL VLDNPHMQQL IQHNPEIGHI LNNPEIMRQT LEFLRNPAMM QEMIRSQDRV LSNLESIPGG YNVLCTMYTD IMDPMLNAVQ EQFGGNPFAT ATTDNATTTT SQPSRMENCD PLPNPWTSTH GGSGSRQGRQ DGDQDAPDIR NRFPNFLGII RLYDYLQQLH ENPQSLGTYL QGTASALSQS QEPPPSVNRV PPSSPSSQEP GSGQPLPEES VAIKGRSSCP AFLRYPTENS TGQGGDQDGA GKSSTGHSTN LPDLVSGGLD SANRVPFAPL SFSPTAAIPG IPEPPWLPSP AYPRSLRPDG MNPAPQLQDE IQPQLPLMH LQAAMANPRA LQALRQIEQG LQVLATEAPR LLLWFMPCLA GTGSVAGGIE SREDPLMSED PLPNPPEVF PALDSAELGF LSPFHLMLQ DLVSTNPQQL QPEAHFQVQL EQLRSMGFLN REANLQALIA TGGDVDAAVE KLRQS
Source	Mammalian cell
Target Names	UBQLN3
Protein Names	Recommended name: Ubiquilin-3
Expression Region	1-655
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	Summary: This gene encodes an ubiquitin-like protein (ubiquilin) that shares high degree of similarity with related products in yeast, rat and frog. Ubiquilins contain a N-terminal ubiquitin-like domain and a C-terminal ubiquitin-associated domain. They physically associate with both proteasomes and ubiquitin ligases, and thus thought to functionally link the ubiquitination machinery to the proteasome to affect in vivo protein degradation. This gene is specifically expressed in the testis, and proposed to regulate cell-cycle progression during



spermatogenesis.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.