



Recombinant Human Rab GTPase-binding effector protein 2 (RABEP2)

Product Code	CSB-EP880980HU-B
Abbreviation	RABEP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H5N1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AAAAPVAAD DDERRRRPGA ALEDSRSQEG ANGEAESGEL SRLRAELAGA LAEMETMKAV AEVSESTKAE AVAAVQRQCQ EEVASLQAIL KDISSYEAQ ITALKQERQQ QQQDCEEKER ELGRLKQLLS RAYPLDSLEK QMEKAHEDSE KLREIVLPM KEIEELKAKL LRAEELIQEI QRRPRHAPSL HGSTELLPLS RDPSPPLEPL EELSGDGGPA AEAFAHNCDD SASISSFSLG GGVGSSSSLP QSRQGLSPEQ EETASLVSTG TLVPEGIYLP PPGYQLVPDT QWEQLQTEGR QLQKDLESVS RERDELQEG LRRSNEDCAKQ MQVLLAQVQN SEQLLRTLQG TVSQAQERVQ LQMAELVTTH KCLHHEVKRL NEENQGLRAE QLPSSAPQGS QQEQGEEESL PSSVPELQQL LCCTRQEARA RLQAQEHGAE RLRIEIVTLR EALEEETVAR ASLEGQLRVQ REETEVLAS LCLSRTEMER VQQEQSKAQL PDLLSEQRAK VLRLQAELET SEQVQRDFVR LSQALQVRLE RIRQAETLEQ VRSIMDEAPL TDVRDIKDT
Source	E.coli
Target Names	RABEP2
Protein Names	Recommended name: Rab GTPase-binding effector protein 2 Alternative name(s): Rabaptin-5beta
Expression Region	2-569
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.