



Recombinant Human UPF0428 protein CXorf56 (CXorf56)

Product Code	CSB-EP880981HU-B
Abbreviation	CXorf56
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H5V9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPKVVSRSVV CSDTRDREEY DDGEKPLHVY YCLCGQMVLV LDCQLEKLPM RPRDRSRVID AAKHAHKFCN TEDEETMYLR RPEGIERQYR KKCAKCGLPL FYQSQPKNAP VTFIVDGAVV KFGQGFGKTN IYTQKQEPPK KVMMTKRTKD MGKFSSVTVS TIDEEEEIE AREVADSYAQ NAKVIEKQLE RKGMSKRRLQ ELAELEAKKA KMKGTLIDNQ FK
Source	E.coli
Target Names	CXorf56
Protein Names	Recommended name: UPF0428 protein CXorf56
Expression Region	1-222
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	While this gene is well-supported by transcript data, no functional information on its protein products is currently available. Three transcript variants encoding different isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.