



Recombinant Human Protein FAM124B (FAM124B)

Product Code	CSB-EP880982HU-B
Abbreviation	FAM124B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H5Z6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDETQGPLAM TVHLLANSGH GSLLQRTLQD LLDICPEVR LFQVSERASP VKYCEKSHSK RSRFPGMSVL LFLHESPGED RLFRVLDSLQ HSPWQCYPTQ DTRGRLCPYF FANQEFYSLD SQLPIWVGRQ VHCSEILRV TLYCSFDNYE DAIRLYEMIL QREATLQKSN FCFVLYASK SFALQLSLKQ LPPGMSVDPK ESSVLQFKVQ EIGQLVPLL NPCMPISSTR WQTQDYDGNK ILLQVQLNPE LGVKNGILGA GMLPLGSRLT SVSAKRTSEP RSQRNQGKRS QGHSLELPEP SGSPTSDRCA GTSWKSPGRS FQVSSPAMGA HLHLSSHLE SGARMKVLNR ENSFQKLEAE TNVDTGLTII NSEPRQTYFG GFPRDLQTSQ PPFCLPASSL GVATSKNNSV LKERVSPLPL AGQRDLGTRK TISECLLHLQ VQGEEKEEDE EEFI
Source	E.coli
Target Names	FAM124B
Protein Names	Recommended name: Protein FAM124B
Expression Region	1-455
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.