



# Recombinant Human Protein FAM124B (FAM124B)

<b>Product Code</b>	CSB-YP880982HU
<b>Abbreviation</b>	FAM124B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H5Z6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDETQGPLAM TVHLLANS GH G SLLQRTL DQ LLD CICEV R LFQV S ERASP VKYCEKSHSK RSRFP GMSVL LFLHESPGED RLFRVLDSLQ HSPWQC YPTQ DTRGR LCPYF FANQEFYSLD SQLPIWGV RQ VHC GSEILRV TLYCSFDNYE DAIRLYEMIL QREATLQKSN FCFFVLYASK SFALQLSLKQ LPPGMSVDPK ESSVLQFKVQ EIGQLVPLL NPCMPISSTR WQTQDYDGNK ILLQVQLNPE LGVKNGILGA GMLPLGSRLT SVSAKRTSEP RSQRNQGKRS QGHSLELPEP SGSPTSDRCA GTSWKSPGRS FQVSSPAMGA HLHLSSHLE SGARMKVLNR ENSFQKLEAE TNVDTGLTII NSEPRQTYFG GFPRDLQTSQ PPFCLPASSL GVATSKNNSV LKERVSPLPL AGQRDLGTRK TISECLLHLQ VQGEEKEEDE EEFFI
<b>Source</b>	Yeast
<b>Target Names</b>	FAM124B
<b>Protein Names</b>	Recommended name: Protein FAM124B
<b>Expression Region</b>	1-455
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.