



Recombinant Human TBC1 domain family member 17 (TBC1D17)

Product Code	CSB-BP881007HU
Abbreviation	TBC1D17
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HA65
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEGAGYRVVF EKGGVYLHTS AKKYQDRDSL IAGVIRVVEK DNDVLLHWAP VEEAGDSTQI LFSKKDSSGG DSCASEEEPT FDPDYEPDWA VISTVRPQLC HSEPTRGAEP SCPQGSWAFS VSLGELKSIR RSKPGLSWAY LVLVTQAGGS LPALHFHRGG TRALLRVLSR YLLASSPQD SRLYLVPDHD SSALSNSFHH LQLFDQDSSN VVSRFLQDPY STTFSSFSRV TNFFRGALQP QPEGAASDLP PPPDEPEPG FEVISCVELG PRPTVERGPP VTEEEWARHV GPEGRLQQVP ELKNRIFSGG LSPSLRREAW KFLGYSWE GTAEHKAHI RKKTDEYFRM KLQWKSVSPE QERRNSLLHG YRSLIERDVS RTDRTNKFYE GPENPGLGLL NDILLTYCMY HFDLGYVQGM SDLLSPILYV IQNEVDAFWC FCGFMELVQG NFEESQETMK RQLGRLLLLL RVLDP LLCDF LDSQDSGSLC FCFRWLLIWF KREFPPDVL RLWEVLWTGL PGPNLHLLVA CAILDMERDT LMLSGFGSNE ILKHINELTM KLSVEDVLR AEALHRQLTA CPELPHNVQE ILGLAPPAEP HSPSPTASPL PLSPTRAPPT PPPSTDTAPQ PDSSLEILPE EEDEGADS
Source	Baculovirus
Target Names	TBC1D17
Protein Names	Recommended name: TBC1 domain family member 17
Expression Region	1-648
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.