



# Recombinant Human TBC1 domain family member 17 (TBC1D17)

<b>Product Code</b>	CSB-EP881007HU-B
<b>Abbreviation</b>	TBC1D17
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9HA65
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEGAGYRVVF EKGGVYLHTS AKKYQDRDSL IAGVIRVVEK DNDVLLHWAP VEEAGDSTQI LFSKKDSSGG DSCASEEEPT FDPDYEPDWA VISTVRPQLC HSEPTRGAEP SCPQGSWAFS VSLGELKSIR RSKPGLSWAY LVLVTQAGGS LPALHFHRGG TRALLRVLSR YLLASSPQD SRLYLVPDHD SSALSNSFHH LQLFDQDSSN VVSRFLQDPY STTFSSFSRV TNFFRGALQP QPEGAASDLP PPPDEPEPG FEVISCVELG PRPTVERGPP VTEEEWARHV GPEGRLQQVP ELKNRIFSGG LSPSLRREAW KFLGYSWE GTAEHKAHI RKKTDEYFRM KLQWKSVSPE QERRNSLLHG YRSLIERDVS RTDRTNKFYE GPENPGLGLL NDILLTYCMY HFDLGYVQGM SDLLSPILYV IQNEVDAFWC FCGFMELVQG NFEESQETMK RQLGRLLLLL RVLDP LLCDF LDSQDSGSLC FCFRWLLIWF KREFPFPDVL RLWEVLWTGL PGPNLHLLVA CAILDMERDT LMLSGFGSNE ILKHINELTM KLSVEDVLTR AEALHRQLTA CPELPHNVQE ILGLAPPAEP HSPSPTASPL PLSPTRAPPT PPPSTDTAPQ PDSSLEILPE EEDEGADS
<b>Source</b>	E.coli
<b>Target Names</b>	TBC1D17
<b>Protein Names</b>	Recommended name: TBC1 domain family member 17
<b>Expression Region</b>	1-648
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.