



Recombinant Human Protein Njmu-R1 (C17orf75)

Product Code	CSB-EP881009HU-B
Abbreviation	C17orf75
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HAS0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLPSLQESMD GDEKELESSE EGGSAEERRL EPPSSSHYCL YSYRGSRLAQ QRGDS EDGSP SGTNAETPSG DDFSLSLADT NLPSEVEPEL RSFIAKRLSR GAVFEG LGNV ASVELKIPGY RVGCYYCLFQ NEKLLPETVT IDSERNPSEY VVCFLGGSEK GLELFRLELD KYIQGLKNNM NCEARGLESH IKSYLSSWFE DVVCPIQRVV LLFQEKL TFL LHAALSYTPV EVKESDEKTK RDINRFLSVA SLQGLIHEGT MTS LCMAMTE EQHKSVIDC SSSQPQFCNA GSNRFCEDWM QAFLNGAKGG NPFLFRQVLE NFKLKAIQDT NNLKRFIRQA EMNHYALFKC YMFLKNCGSG DILLKIVKVE HEEMPEAKNV IAVLEEFMKE ALDQSF
Source	E.coli
Target Names	C17orf75
Protein Names	Recommended name: Protein Njmu-R1
Expression Region	1-396
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.