



Recombinant Human E3 ubiquitin-protein ligase pellino homolog 2 (PELI2)

Product Code	CSB-YP881010HU
Abbreviation	PELI2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HAT8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MFSPGQEEHC APNKEPVKYG ELVVLYNGA LPNGDRGRRK SRFALYKRPK ANGVKPSTVH VISTPQASKA ISCKGQHSIS YTLARNQTVV VEYTHDKDTD MFQVGRSTES PIDFVVTDTI SGSQNTDEAQ ITQSTISRFA CRIVCDRNEP YTARIFAAGF DSSKNIFLGE KAAKWKNPDG HMDGLTTNGV LVMHPRGGFT EESQPGVWRE ISVCGDVYTL RETRSAQQRG KLVESETNVL QDGLIDLGC ATLLWRTADG LFHTPTQKHI EALRQEINAA RPQCPVGLNT LAFPSINRKE VVEEKQPMWAY LSCGHVHGYH NWGHRSDTEA NERECPCMRT VGPYVPLWLG CEAGFYVDAG PPTHAFTPCG HVCSEKSAKY WSQIPLPHGT HAFHAACPFC ATQLVGEQNC IKLIFQGPID
Source	Yeast
Target Names	PELI2
Protein Names	Recommended name: E3 ubiquitin-protein ligase pellino homolog 2 Short name= Pellino-2 EC= 6.3.2.-
Expression Region	1-420
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.