



# Recombinant Human Netrin-4 (NTN4)

<b>Product Code</b>	CSB-MP881014HU
<b>Abbreviation</b>	NTN4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9HB63
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	GL SGVAGVSSRC EKACNPRMGN LALGRKLWAD TTCGQNATEL YCFYSENTDL TCRQPKCDKC NAAYPHLAHL PSAMADSSFR FPRTWWQSAE DVHREKIQLD LEAEFYFTHL IVMFKSPRPA AMVLDRSQDF GKTWKPKYKYF ATNCSATFGL EDDVVKKGAI CTSKYSSPFP CTGGEVIFKA LSPPYDTENP YSAKVQEQLK ITNLRVQLLK RQSCPCQRND LNEEPQHFTH YAIYDFIVKG SCFCNGHADQ CIPVHGFRPV KAPGTFHMHVH GKCMCKHNTA GSHCQHCAPL YNDRPWEAAD GKTGAPNECR TCKCNGHADT CHFDVNVWEA SGNRSGGVCD DCQHNTGQY CQRCKPGFYR DLRRPFSAPD ACKPCSCHPV GSAVLPANSV TFCDPSNGDC PCKPGVAGRR CDRCMVGYWG FGDYGCRPCD CAGSCDPITG DCISSHTDID WYHEVPDFRP VHNKSEPAWE WEDAQGFSAL LHSGKCECKE QTLGNAKAFK GMKYSYVLKI KILSAHDKGT HVEVNVKIKK VLKSTKLKIF RGKRTLYPES WDRGCTCPI LNPGLYLVLA GHEDIRTGKL IVNMKSFVQH WKPSLGRKVM DILKRECK
<b>Source</b>	Mammalian cell
<b>Target Names</b>	NTN4
<b>Protein Names</b>	Recommended name: Netrin-4 Alternative name(s): Beta-netrin Hepar-derived netrin-like protein
<b>Expression Region</b>	19-628
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.