



Recombinant Human Kinetochore protein Spc25 (SPC25)

Product Code	CSB-BP881019HU
Abbreviation	SPC25
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HBM1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVEDELALFD KSINEFWNKF KSTDTSCQMA GLRDTYKDSI KAFAEKLSVK LKEEERMVEM FLEYQNQISR QNKLIQEKKD NLLKLAIEVK GKKQELEVLT ANIQDLKEEY SRKKETISTA NKANAERLKR LQKSADLYKD RLGLEIRKIY GEKLQFIFTN IDPKNPESPF MFSLHLNEAR DYEVS DSAPH LEGLAEFQEN VRKTNNFSAF LANVRKAFTA TVYN
Source	Baculovirus
Target Names	SPC25
Protein Names	Recommended name: Kinetochore protein Spc25 Short name= hSpc25
Expression Region	1-224
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a protein that may be involved in kinetochore-microtubule interaction and spindle checkpoint activity.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.