



Recombinant Danio rerio CUGBP Elav-like family member 1 (celf1), partial

Product Code	CSB-EP881258DIL-B
Abbreviation	celf1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9IBD0
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MNGSLDHPDQ PDIDSIKMFV GQIPRTWSED QLRELFEPYG AVYEINVLRD RSQNPPQSKG CCFVTTYTRK SALEAQNALH NMKILPGMHM PIQMKPADSE KNAVEDRKL FVGMISKKCN ENDIRLMFSP YGQIEECRIL RGPDGLSRGC AFVTFTARQM AQSIAKSMHQ SQTMEGCSSP IVVKFADTQK DKEQKRIAQQ LQQQMQLNA ASMWGNLTGL NSLGPQYLAL YLQLLQQSAS SGNALNNLHP MSGLNAMQNL AALAAAASAT QATPTGSSAL TTSSSPLSVL TSSGTPSGQP AQSAWDAYKA GSSPTSSTSS SVNPMASLGA LQSLAAGAGA GLNMSSLASM AALNGGLGSG GLSNGSGSTM EALTQAAYSG IQQYAAAALP SLYSQSLLSQ QNVSAAGSQK EGPEGANLFI YHLPQEFGDQ DLLQMFMPFG NVISAKVFID KQTNLKCFG FVSYDNPVSS QAAIQSMNGF QIGMKRLKVQ LKRSKNDSP Y
Source	E.coli
Target Names	celf1
Protein Names	Recommended name: CUGBP Elav-like family member 1 Short name= CELF-1 Alternative name(s): Bruno-like protein 2 CUG triplet repeat RNA-binding protein 1 Short name= CUG-BP1 CUG-BP- and ETR-3-like factor 1 EDEN-BP/Bruno-lik
Expression Region	1-501
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.