



# Recombinant Mouse Cullin-3 (Cul3)

<b>Product Code</b>	CSB-BP881350MO
<b>Abbreviation</b>	Cul3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9JLV5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	SNLSKGTGS RKDTKMRIRA FPMTMDEKYV NSIWDLLKNA IQEIQRKNNS GLSFEELYRN AYTMLHKG EKLYTGLREV VTEHLINKVR EDVLNSLNNN FLQTLNQAWN DHQTAMVMIR DILMYMDRVY VQQNNVENVY NLGLIIFRDQ VVRYGCIRDH LRQTLDMIA RERKGEVVDR GAIRNACQML MILGLEGRSV YEEDFEAPFL EMSAEFFQME SQKFLAENSA SVYIKKVEAR INEEIERVMH CLDKSTEEPI VKVVERELIS KHMKTIVEME NSGLVHMLKN GKTEDLACMY KLFSRVPNGL KTMCECMSCY LREQGKALVS EEGEGKNPVD YIQGLLDLKS RFDRFLQESF NNDRLFQTI AGDFEYFLNL NSRSPEYLSL FIDDKLKKGV KGLTEQEVET ILDKAMVLFV FMQEKDVFER YYKQHLARRL LTNKSVSDDS EKNMISKLKT ECGCQFTSKL EGMFRDMSIS NTTMDEFQRH LQATGVSLGG VDLTVRVLTT GYWPTQSATP KCNIPPAPRH AFEIFRRFYL AKHSGRQLTL QHMGASADLN ATFYGPVKKE DGSEVGVGGA QVTGSNTRKH ILQVSTFQMT ILMLFNNREK YTFEEIQQET DIPERELVRA LQSLACGKPT QRVLTKEPKS KEIESGHIFT VNDQFTSKLH RVKIQTVAQK QGESDPERKE TRQKVDDDRK HEIEAAIVRI MKSRKKMQHN VLVAEVTQQL KARFLPSPVV IKKRIEGLIE REYLARTPED RKVYTYVA
<b>Source</b>	Baculovirus
<b>Target Names</b>	Cul3
<b>Protein Names</b>	Recommended name: Cullin-3 Short name= CUL-3
<b>Expression Region</b>	2-768
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.