



Recombinant *Synechococcus elongatus* Quinolinate synthase A (nadA)

Product Code	CSB-YP881687FPY
Abbreviation	nadA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9L790
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	>85% (SDS-PAGE)
Sequence	MFLAADRPTT DLPADLPAAI AALKQELNAV ILAHYYQEAA IQDVADYIGD SLGLSRQAAA ADADVIVFAG VHFMAETAKI LNPQRQVLLP DLAAGCSLAD SCPPEAFAAF KAAHPNHIVI SYINCTAEIK ALSDICTSS NAVKIVQQIP VDQPIIFAPD RNLGRYVMQQ TGRDLVLWDG SCIVHETFSE QRLLLELQARH PDAEIIAHPE CETPVLQAR FIGSTTALLN YSLNSPSREF IVVTEPGIIH QMQQAPEKT FIPAPPQDTT CACNECPFMR LNTLEKLYLC MRDRRPEIQI PEETRLAALR PIERMLAMSA
Source	Yeast
Target Names	nadA
Protein Names	Recommended name: Quinolinate synthase A EC= 2.5.1.72
Expression Region	1-320
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.