



# Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At5g14820, mitochondrial (At5g14820)

<b>Product Code</b>	CSB-EP881710DOA
<b>Abbreviation</b>	At5g14820
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9LEQ7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SD GCDEEVESEC DNDEETGVSC VESSTNPEEV ERVCKVIDEL FALDRNMEAV LDEMKLDSLH DLIVEVLERF RHARKPAFRF FCWAAERQGF AHDSRTYNSM MSILAKTRQF ETMVSVLEEM GTKGLLTMET FTIAMKAFAA AKERKKA VGI FELMKKYKFK IGVETINCLL DSLGRAKLGK EAQVLFDKLLK ERFTPNMPTY TVLLNGWCRV RNLIEAARIW NDMIDHGLKP DIVAHNVMLE GLLRSMKKSD AIKLFHVMKS KGPCPNVRSY TIMIRDFCKQ SSMETAIEYF DDMVDSGLQP DAAVYTCLIT GFGTQKKLDT VYELLKEMQE KGHPPDGKTY NALIKLMANQ KMPEHGTRIY NKMIQNEIEP SIHTFNMIMK SYFVARNYEM GRAVWDEMIK KGICPDDNSY TVLIRGLISE GKSREACRYL EEMLDKGMKT PLIDYNKFAA DFHRGGQPEI FEELAQRAKF SGKFAAAEIF ARWAQMTRRR CKQRFMED
<b>Source</b>	E.coli
<b>Target Names</b>	At5g14820
<b>Protein Names</b>	Recommended name: Pentatricopeptide repeat-containing protein At5g14820, mitochondrial
<b>Expression Region</b>	99-598
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . The shelf life of lyophilized form is 12 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ .