



Recombinant Arabidopsis thaliana F-box protein PP2-A15 (PP2A15)

Product Code	CSB-EP881716DOA-B
Abbreviation	PP2A15
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LF92
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGSSLSNLND GTNGLAMGPG LGDIPESCVV CVMYLTPE ICNLAGLNRS FRGAASSDSV WEKKLPENYQ DLLDLLPPER YHLSKKDIF AVLSRPIPFD DDNKEVWIDR VTGRVCMAS ARGMSITGIE DRRYWNWIPT EESRFHVVAY LQQIWWFEVD GTVRFHLPPG VYSLSFRIHL GRFTKRLGRR VCHFELTHGW DLKPVRFSLS TSDGQEASCE YYLDDVERNE ALGKHKRGYW IEYRVGEFIV NGSEPSTEIQ WSMKQIDCTH SKGGLCVDSV FINPIGDVKE HKRKAVLKEA
Source	E.coli
Target Names	PP2A15
Protein Names	Recommended name: F-box protein PP2-A15 Alternative name(s): Protein PHLOEM PROTEIN 2-LIKE A15 Short name= AtPP2-A15
Expression Region	1-300
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.