



# Recombinant Arabidopsis thaliana F-box/kelch-repeat protein At3g13680 (At3g13680)

<b>Product Code</b>	CSB-EP881735DOA
<b>Abbreviation</b>	At3g13680
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9LID1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTTMGDLPGD LVEEILSRVP LTSLRAIRST CQKWNSLSKS QICGRKATAA ENKFLGFMMK DSRVCSMKFD LQGIRNDDGE LVEPSIKQVS KLDQIEVSQV FHCDGLVLCI IKDNTGLLVW NPYLGQTRWI HPRNNYQIED RYALGYDNNR NHKILRIFDL YPSSNRVFGY EYDFSSNSW KVL DVIPEWD IQSHFRGASL KGNAYFPAQK KETVGGIKTI NIEDVLLCFD FTRERFGPPL PLPFHSYNAD FFVSLACVRE EQLAVLYQRW GAFETIEICV TNKIDPNTVS WSKFLITSTG FPVDTFSGSF FIDEEKKVA VFDLDRYKPT ETCRYQIVYI IGQDEYFKSV NMGVAPNIAK NSYQAYCVPL VCSSSYVPSL VQLQINKPGK RKESN
<b>Source</b>	E.coli
<b>Target Names</b>	At3g13680
<b>Protein Names</b>	Recommended name: F-box/kelch-repeat protein At3g13680
<b>Expression Region</b>	1-395
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.