



Recombinant Arabidopsis thaliana PHD finger protein ING1 (ING1)

Product Code	CSB-BP881738DOA
Abbreviation	ING1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LIQ6
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MSFAEEFEAN LVSLAHVLQK KYALLRDLK SLQENQRQNE QRCEKEIEDI RRGRAGNITP NTS�TKFSEE ALDEQKHSVR IADEKVTLAM QAYDLVDMHV QQLDQYMKKS DEVIRKEKEA AAATLELENN GKAGNAGEGG RGGRKKTRLA TAASTAAST GMTSSNMDLD LPVDPNEPTY CİCNQVSFGE MVACDNNACK IEWFHFGCVG LKEQPKGKWKY CPECATVKKS RKGR
Source	Baculovirus
Target Names	ING1
Protein Names	Recommended name: PHD finger protein ING1 Alternative name(s): Protein INHIBITOR OF GROWTH 1 Short name= Protein AtING1
Expression Region	1-234
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.