



Recombinant Arabidopsis thaliana Putative F-box protein At3g20030 (At3g20030)

Product Code	CSB-EP881744DOA
Abbreviation	At3g20030
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LJZ8
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MTMMSDLSQD LLEEILSRVP RTSLGAVRST CKRWNTLFKD RILCKAEETR DQFRFIMKKY KLCSMRFDLN GTLNEDGGTE FVDPSIKELG HFFNQVKISK VFQCDGLLLC VTKEDNIRLV VWNPYLGQIR WIESGNTNYR LFDYAIGYD NNRKHKILRF LEYYNFMKH RFLEYEYDF SSNSWRVLDI APRWEIESYQ RGASLKGNTY FIAKEKIEYE EDGEFPEPAD NLLCFDFTTE SFGQFLPLPF QHLYDVGAL SSLGDEKLAA LFCQGD TDLS EVEIWVTTLT ETNTVSWNPF LKVDMEPHYG RSFMDYGG SFFIDEEKKL AVVFQFDESG MTRYEDATYI IGENGYVKKV RLGEAPANQG GYCFPSVCFS SYVPSLVQIN QIVGFKKERE EE
Source	E.coli
Target Names	At3g20030
Protein Names	Recommended name: Putative F-box protein At3g20030
Expression Region	1-402
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.