



# Recombinant Arabidopsis thaliana Putative F-box protein At1g12855 (At1g12855)

<b>Product Code</b>	CSB-EP881774DOA-B
<b>Abbreviation</b>	At1g12855
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9LPW4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MESREDSFIS KEKKSTMKKE KQAIASQRNR RRVIKNRGNG KRLIASLSQR KRRRIPRGRG NEKAVFAPSS LPNDVVEEIF LRLPVKAIQ LKSLSKQWRS TIESRSFEER HLKIVERSRV DFPQVMVMSE EYSLKGSKGN QPRPDTDIGF STICLESASI LSSTLITFPQ GFQHRIYASE SCDGLFCIHS LKTQAIYVVN PATRWFQRQLP PARFQILMQK LYPTQDTWID IKPVVCYTAF VKANDYKLVW LYNSDASNPV LGVTKCEVFD FRANAWRYLT CTPSYRIFPD QVPAATNGSI YWFTEPYNGE IKVVALDIHT ETRVLPKIN PAIASSDPDH IDMCTLDNGL CMSKRESDTL VQEIWRLKSS EDSWEKVYTI DLLSCSSSSL SEFRDGFNWT RKDLVEPSTP VAICKDKKIL LSHRYARNMI KYDPQIKSIS LIYQPPLCRR YASYFQSLIS HI
<b>Source</b>	E.coli
<b>Target Names</b>	At1g12855
<b>Protein Names</b>	Recommended name: Putative F-box protein At1g12855
<b>Expression Region</b>	1-462
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.