



Recombinant Arabidopsis thaliana Putative HVA22-like protein g (HVA22G)

Product Code	CSB-MP881781DOA
Abbreviation	HVA22G
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LR09
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MIGSFLTRGL LMFVGYAYPA YECFKTVELN KPEIQQLQFW CQYWIIVAAL TIFERIGDAL VSWLPMYSEA KLAFFIYLWF PKTKGTTYVY DSFFRPYIAK HENEIDRNLV KVKTRAKDMA MIYLQKAINQ GQTKFFEILQ YITEQSTPKS KAEEKKETTI PKLDDPILKV KENEVTK
Source	Mammalian cell
Target Names	HVA22G
Protein Names	Recommended name: Putative HVA22-like protein g Short name= AtHVA22g
Expression Region	1-177
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.