



Recombinant Arabidopsis thaliana Probable E3 ubiquitin protein ligase DRIPH (At3g23060)

Product Code	CSB-BP881793DOA
Abbreviation	At3g23060
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LS86
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MLTKVLSKEV KPCLACPICT NPFKDATTIS ECLHTFCRSC IRNKFINERV NACPVCNVNL GVFPLEKLR S DCTWQDLK LK IYRAMMESLK KAGPKTVAAS VKSSKKKRKS RTSLRVSSSR VSSSPDTPLE PANVVVEPPN VVVEEKHRET VLALQSTRKP IITFQKRGRK ASLPKKIDSK PEPELPPKEP KIKNLFDLNN EPEDNGLDEA EGSTFQEVVP KEKDLCKPIF SLSVTLNIND TPPDIVEPEI SSDDDTEESV EPIQNKCVVN RETKEVPVQV NQNSLLISSD RDREDNSGQK LKTNGAATSR SRKKKGKPV EKSYSLRPRI GRRTVNPAAG TTTPEAPVSV EEEMKVEEGR NNNPVWFSK PSKTQNIEM LPPITACCIR VKDSNMTVSY LKKYLMVKLG LESEDQVEIW LRDEPVCSSL TLHNLVDWWW QTTP LPERQS AMVGSSAAEF IMDLYYSFKS DASDSGSASE
Source	Baculovirus
Target Names	At3g23060
Protein Names	Recommended name: Probable E3 ubiquitin protein ligase DRIPH EC= 6.3.2.- Alternative name(s): DREB2A-interacting protein homolog
Expression Region	1-480
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.