



Recombinant Arabidopsis thaliana Probable E3 ubiquitin protein ligase DRIPH (At3g23060)

| | |
|--------------------------|--|
| Product Code | CSB-EP881793DOA |
| Abbreviation | At3g23060 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9LS86 |
| Product Type | Recombinant Protein |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MLTKVLSKEV KPCLACPICT NPFKDATTIS ECLHTFCRSC IRNKFINERV NACPVCNVNL GVFPLEKLR S DCTWQDLK LK IYRAMMESLK KAGPKTVAAS VKSSKKKRKS RTSLRVSSSR VSSSPDTPLE PANVVVEPPN VVVEEKHRET VLALQSTRKP IITFQKRGRK ASLPKKIDSK PEPELPPKEP KIKNLFDLNN EPEDNGLDEA EGSTFQEVVP KEKDLCKPIF SLSVTLNIND TPPDIVEPEI SSDDDTEESV EPIQNKCVVN RETKEVPVQV NQNSLLISSD RDREDNSGQK LKTNGAATSR SRKKKGKPV EKSYSLRPRI GRRTVNPAAG TTTPEAPVSV EEEMKVEEGR NNNPVWFSK PSKTQNIEM LPPITACCIR VKDSNMTVSY LKKYLMVKLG LESEDQVEIW LRDEPVCSSL TLHNLVDWWV QTTP LPERQS AMVGSSAAEF IMDLYYSFKS DASDSGSASE |
| Source | E.coli |
| Target Names | At3g23060 |
| Protein Names | Recommended name: Probable E3 ubiquitin protein ligase DRIPH EC= 6.3.2.- Alternative name(s): DREB2A-interacting protein homolog |
| Expression Region | 1-480 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.