



Recombinant Arabidopsis thaliana B3 domain-containing protein At3g25182 (At3g25182)

Product Code	CSB-MP881797DOA
Abbreviation	At3g25182
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LSF7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MTSMYYNHHDH LTASGKNMWS NFYVLVDTAV MLYEEEEKQRR KIVSEEEEEET QKRIFCLFPR KTRSSLVKRQ QKLSGFLTSA SSSLIDLNQF PSDSEIEDPL NHHQLLSSSC FIAADFETLQ NPSSEPCSSL VLFDCCKTAGS EKTETKDPPN PNFPCPLSLC LTENKNRKRR AVEQRKRTGG VKKAKVAPFS QTARETPEWL VNVMRDMKEA KDAKLIFEKT LFTVDVNPTQ NRLSMPFNLL LRNDFLTSVE SRIIDKDIKN DKKIGVGAIL VDQRCKKWGV MLKRWEMKKE SGKGSWNYNL ICGWNDIVEA NGLKEGDNIS LWSFRSCGIL CFAMEQRSPS LALCLC
Source	Mammalian cell
Target Names	At3g25182
Protein Names	Recommended name: B3 domain-containing protein At3g25182
Expression Region	1-346
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.