



Recombinant Arabidopsis thaliana Putative F-box/LRR-repeat protein At3g28410 (At3g28410)

Product Code	CSB-MP881798DOA
Abbreviation	At3g28410
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LSJ3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MAKRREKRGR RRQHRSHRKI QRIIDGADFI NYMPDDILHH ILSFIPTDLA MRTSVLSRRW RHVWCETPCL DIKLKHGETN QTLTSYTAPI ITSFKLVMDL NDNTVPQVDS WIEFALSRNV QNLSVFVRDF TYTKTYRFPD IFYISSSLKQ LDVTLDDFFDM IPTCAVSWKS LRNLTLRFCQ IPDESMHNIL SGCPILESLT LDTCRLLERL DLSKSPNLRR LDINRQYRRT GPIAIVAPHI YYLRLTYSST PSTIVDVSSL SEANLNII SD RLLSPLTADR YQTMAL EMLS KFHNVKRLTV GETILQILSL AELRGVPFPT LKVQTLTVKT EFVRSVIPGI SROLLQNSPGL KKLTLHTLQL SHDIMEMHPL KGLYPDQCWR STCEVFPTSK EIYKMLGCND ATSKLVASF M NLVLRNAKTL ERMVVWLGGI YFNGDAPWFE EELFDMVETL SHNNNVSILL KQSNC
Source	Mammalian cell
Target Names	At3g28410
Protein Names	Recommended name: Putative F-box/LRR-repeat protein At3g28410
Expression Region	1-465
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.