



Recombinant Arabidopsis thaliana Putative F-box protein At3g17490 (At3g17490)

Product Code	CSB-EP881827DOA
Abbreviation	At3g17490
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LUP8
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MMMPHLSIDL VEEILSRVPA ISLKRLRYTC KQWNALFNDQ RFSKKHRDKA PKTYLGLTLK DFRIYSMSSN LHGLLHNNNI DLLMEFKGKL SSLNDLNDFE ISQIYPCDGL ILCSTKRNTR LVVWNPCTGQ TRWIKRRNRR MCDTFAFGYD NSKSSCLNNY KILRVCEKIK GQQFEYEIFE FSSNSWRVLD VNPNCIIEGR SVSVKGNYSW FATITKTHYF IRRFDFSSET FQKLPLPFHI FDYNDLRALS AFREEQLSVL HQSFDTEKMD IWVTNKIDET TDWSWSKFFT VRLINRLDYP ISMMMKSPLS FFIDEKKNII LCYDKHRENT YKSLVLIVGK DKVYREFYFP ESYELGRTYL CNYVPSLVQI KQSGLISTRK RKRKDFSS
Source	E.coli
Target Names	At3g17490
Protein Names	Recommended name: Putative F-box protein At3g17490
Expression Region	1-388
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.