



Recombinant Arabidopsis thaliana Protein yippee-like At3g55890 (At3g55890)

Product Code	CSB-BP881859DOA
Abbreviation	At3g55890
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LY56
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGRVFMVDLE GNIYICKLCK THLSTDQDIM SKSFQCKNGR AYLFNNVVNV SVGEKEDRMM ITGLHNVVDI FCVGC GSNVG WKYEFAHEKS QKYKEGKSVL ELYKISGPHD SNDLVSDGDD A
Source	Baculovirus
Target Names	At3g55890
Protein Names	Recommended name: Protein yippee-like At3g55890
Expression Region	1-121
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.