



Recombinant Petunia hybrida Chalcone--flavonone isomerase C (CHI3)

Product Code	CSB-MP881919EUD
Abbreviation	CHI3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9M5B3
Product Type	Recombinant Protein
Immunogen Species	Petunia hybrida (Petunia)
Purity	>85% (SDS-PAGE)
Sequence	MSPPVSVTKM QVENYAFAPT VNPAGSTNTL FLAGAGHRGL EIEGKFVKFT AIGVYLEESA IPFLAEKWKG KTPQELTDSV EFFRDVVTGP FEKFTRVTTI LPLTGKQYSE KVAENCVAHW KGIGTYTDDE GRAIEKFLDV FRSETFPPGA SIMFTQSPLG LLTISFAKDD SVTGTANAVI ENKQLSEAVL ESIIGKHGVS PAAKCNVAER VAELLKKSya EEASVFGKPE TEKSTIPVIG V
Source	Mammalian cell
Target Names	CHI3
Protein Names	Recommended name: Chalcone--flavonone isomerase C Short name= CHI-C Short name= Chalcone isomerase C EC= 5.5.1.6
Expression Region	1-241
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.