



Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At1g05750, chloroplastic (PDE247)

Product Code	CSB-BP881941DOA
Abbreviation	PDE247
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9MA50
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	AAKEFS DMTLAGVEPN HITFIALLSG CGDFTSGSEA LGDLLHGYAC KLGLDRNHVM VGTAIIGMYS KRGRFKKARL VFDYMEDKNS VTWNTMIDGY MRSGQVDNAA KMFDKMPERD LISWTAMING FVKKGYQEEA LLWFREMQIS GVKPDYVAII AALNACTNLG ALSFGLWVHR YVLSQDFKNN VRVSNLIDL YCRCGCVEFA RQVFYNMEKR TVVSWNSVIV GFAANGNAHE SLVYFRKMQE KGFKPDAVTF TGALTACSHV GLVEEGLRYF QIMKCDYRIS PRIEHYGCLV DLYSRAGRLE DALKLVQSMP MKPNEVVIGS LLAACSNHGN NIVLAERLMK HLTDLNVKSH SNYVILSNMY AADGKWE GAS KMRRKMKGLG LKKQPGFSSI EIDDCMHVFM AGDNAHVETT YIREVLELIS SDLRLQGCVV ETLAGDLLNA
Source	Baculovirus
Target Names	PDE247
Protein Names	Recommended name: Pentatricopeptide repeat-containing protein At1g05750, chloroplastic Alternative name(s): Protein PIGMENT DEFECTIVE 247
Expression Region	55-500
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.