



Recombinant Bovine Voltage-dependent L-type calcium channel subunit beta-3 (CACNB3)

| | |
|--------------------------|---|
| Product Code | CSB-YP882003BO |
| Abbreviation | CACNB3 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9MZL3 |
| Product Type | Recombinant Protein |
| Immunogen Species | Bos taurus (Bovine) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MYDDSYVPGF EDSEAGSADS YTSRPSLDSD VSLEEDRESA RREVESQAQQ QLERAKHKPV AFAVRTNVSY CGVLDEECPV QGSGVNFEAK DFLHIKEKYS NDWWIGRLVK EGGDIAFIPS PQRLESIRLK QEQAARRSGN PSSLSDIGSR RSPPPSLAKQ KQKQAEHIPP YDVVPSMRPV VLVGPSLKG Y EVTDMMQKAL FDLKHFRFDG RISITRVTAD LSLAKRSVLN NPGKRTHIER SSARSSIAEV QSEIERIFEL AKSLQLVVDL ADTINHPAQL AKTSLAPIV FVKVSSPKVL QRLIRSRGKS QMKHLTVQMM AYDKLVQCPP ESFDVILDEN QLEDACEHLA EYLEVYWRAT HHPAPGPGLL GPPSAIPGLQ NQQLLGERGE EHSPLERDSL MPSDEASESS RQAWTGSSQR SSRHLEEDYA DAYHDLYQPH RQHTSGLPSA NGHDPQDRLL AQDSEHNHNE RNWQRNRPWP KDSY |
| Source | Yeast |
| Target Names | CACNB3 |
| Protein Names | Recommended name: Voltage-dependent L-type calcium channel subunit beta-3 Short name= CAB3 Alternative name(s): Calcium channel voltage-dependent subunit beta 3 |
| Expression Region | 1-484 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.