



# Recombinant Uncharacterized protein Y39B6A.33 (Y39B6A.33)

<b>Product Code</b>	CSB-EP882041CXY-B
<b>Abbreviation</b>	Y39B6A.33
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NEU5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVAGKRTGAA KGSRHNKKYW RKGNTNIDDIE DSIHIKSRQA ATGGVISEMK DEDLFIVDRT ATANKPVVPK LTKKQQAAL KITKNITQEH VTLPKPSTTS KILKKPAKLP RGNAILALKK GPKAAAPAAK KKNFDVWTTD LTPKIPKSKL ENQEAAEHFL KVVKKKQPKT PGKSITSLLP AVQIAEGGAS YNPESAHEYQE YVAKIAGEEQ KLIDHEAKIK AGIEPQWEKV TTEHERFLEM AEGLRHPKY GKDDEEEEEEA GNSEKSMKTG GEAEPKSQRV ECDRMTKEQK KKKAKAQKLD KEEKRRLEEK AKEQDSHNVY RTKQLHKELD EEEKQRHEES EVRKKEKLN KLTKRQQLGK GKFVDAEDPF LLQEELTGNL RQLKPQGHVL DDRMKSLQRR NMLPIGGDKE KRRIKNRLKS KVV EKRS AKN IVKGS RVI
<b>Source</b>	E.coli
<b>Target Names</b>	Y39B6A.33
<b>Protein Names</b>	Recommended name: Uncharacterized protein Y39B6A.33
<b>Expression Region</b>	1-438
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.