



Recombinant *Drosophila yakuba* Trithorax group protein osa (osa)

Product Code	CSB-EP882045DMR-B
Abbreviation	osa
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NGB4
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila yakuba</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	YLAAADSAMA RTVALQSPCI SYLVAFIEQA EQTALGVANQ HGINYLRENPD DSMGTSLDML RRAAGTLLHL AKHPDNRSLF MQQEQRLLGL VMSHILDQQV ALIISRVLYQ VSRGTGPIHS VEFRLQQRQ QQQLRPASAE KQAASAGGSA PVKAEAASTE TSSTEAKPAP AATTAVVNDE NSNSSQQLPP AATFNDVSN STNSNSCGTV SSNQTNXXX NSSHSSSAIS SQSAITVTAP SAPATGATSA XXXAITSDDQ QVSKVAAAAA AAAALSNASA AAAAAAAAAA ASVGPPTSS VSAGAAVAQP AAPPTNAGT TTAV
Source	E.coli
Target Names	osa
Protein Names	Recommended name: Trithorax group protein osa Alternative name(s): Protein eyelid
Expression Region	1-324
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.