



Recombinant Human Vacuolar protein sorting-associated protein VTA1 homolog (VTA1)

| | |
|--------------------------|---|
| Product Code | CSB-YP882062HU |
| Abbreviation | VTA1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9NP79 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | AALAPLPPL PAQFKSIQHH LRTAQEHDKR DPVVAYYCRL YAMQTGMKID SKTPECRKFL SKLMDQLEAL KKQLGDNEAI TQEIVGCAHL ENYALKMFLY ADNEDRAGR F HKNMKSFYT ASLLIDVITV FGELTDENVK HRKYARWKAT YIHNCLKNGE TPQAGPVGIE EDNDIEENED AGAASLPTQP TQPSSSSTYD PSNMPSGNYT GIQIPPGAHA PANTPAEVPH STGVASNTIQ PTPQTIPAI PALFNTISQG DVRLTPEDFA RAQKYCKYAG SALQYEDVST AVQNLQKALK LLTTGRE |
| Source | Yeast |
| Target Names | VTA1 |
| Protein Names | Recommended name: Vacuolar protein sorting-associated protein VTA1 homolog Alternative name(s): Dopamine-responsive gene 1 protein Short name= DRG-1 LYST-interacting protein 5 Short name= LIP5 SKD1-binding protein 1 Short |
| Expression Region | 2-307 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.