



Recombinant Human Ras-related protein Rab-9B (RAB9B)

Product Code	CSB-YP882063HU
Abbreviation	RAB9B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NP90
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSGKSLLLKV ILLGDGGVGK SSLMNRVYVTN KFDSQAFHTI GVEFLNRDLE VDGRFVTLQI WDTAGQERFK SLRTPFYRGA DCCLLTFSVD DRQSFENLGN WQKEFIYYAD VKDPEHFV VLGNKVDKED RQVTTEEAQT WCMENGDYPY LETSAKDDTN VTVAFEEAVR QVLAVEEQLE HCMLGHTIDL NSGSKAGSSC C
Source	Yeast
Target Names	RAB9B
Protein Names	Recommended name: Ras-related protein Rab-9B Alternative name(s): Rab-9-like protein Short name= Rab-9L
Expression Region	1-201
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.