



# Recombinant Human PR domain-containing protein 11 (PRDM11)

<b>Product Code</b>	CSB-EP882077HU-B
<b>Abbreviation</b>	PRDM11
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NQV5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLKMAEPIAS LMIVECRACL RCSPLFLYQR EKDRMTENMK ECLAQTNAAV GDMVTVVKTE VCSPLRDQEY GQPCSRRPDS SAMEVEPKKL KGKRDIVPK SFQQVDFWFC ESCQEYFVDE CPNHGPPVFV SDTPVPVGIP DRAALTIPQG MEVVKDTSGE SDVRCVNEVI PKGHIFGPYE GQISTQDKSA GFFSWLIVDK NNRYKSIDGS DETKANWMRY VVISREEREQ NLLAFQHSEY IYFRACRDIR PGEWLRVWYS EDYMKRLHSM SQETIHRNLA RGEKRLQREK SEQVLDNPED LRGPIHLSVL RQGKSPYKRG FDEGDVHPQA KKKKIDLIFK DVLEASLESA KVEAHQLALS TSLVIRKVPK YQDDAYSQCA TTMTHGVQNI GQTQGEQDWK VPQGVSKPEG QLEDEEEEPS SFKADSPAEA SLASDPHELP TTSFCPCNIR LKKKVRELQA ELDMLKSGKL PEPPVLPPQV LELPEFSDPA GKLVWMRLLS EGRVRSGLCG G
<b>Source</b>	E.coli
<b>Target Names</b>	PRDM11
<b>Protein Names</b>	Recommended name: PR domain-containing protein 11
<b>Expression Region</b>	1-511
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.