



Recombinant Human Vacuolar protein sorting-associated protein 45 (VPS45)

Product Code	CSB-BP882092HU
Abbreviation	VPS45
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NRW7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MNVVFAVKQY ISKMIEDSGP GMKVLLMDKE TTGIVSMVYT QSEILQKEVY LFERIDSQNR EIMKHLKAIC FLRPTKENVD YIIQELRRPK YTIYFIYFSN VISKSDVKSL AEADEQEVVA EVQEFYGDYI AVNPHLFLSLN ILGCCQGRNW DPAQLSRTTQ GLTALLLSLK KCPMIRYQLS SEAAKRLAEC VKQVITKEYE LFEFRRTEVP PLLLILDRCD DAITPLLQW TYQAMVHELL GINNNRIDLS RVPGISKDLR EVVLSAENDE FYANNMYLNF AEIGSNIKNL MEDFQKKKPK EQQKLESIAD MKAFVENYPQ FKKMSGTVSK HVTVVGELSR LVSERNLLEV SEVEQELACQ NDHSSALQNI KRLLQNPKVT EFDAARLVML YALHYERHSS NSLPGLMMDL RNKGVSEKYR KLVS AVVEYG GKRVRGSDLF SPKDAVAITK QFLKGLKGVE NVYTQHQPFL HETLDHLIKG RLKENLYPYL GPSTLRDRPQ DIIVFVIGGA TYEEALTVYN LNRTTPGVRI VLG GTTVHNT KSFLEEVLAS GLHSRSKESQ QVTSRSASRR
Source	Baculovirus
Target Names	VPS45
Protein Names	Recommended name: Vacuolar protein sorting-associated protein 45 Short name= h-VPS45 Short name= hVps45
Expression Region	1-570
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	Vesicle mediated protein sorting plays an important role in segregation of intracellular molecules into distinct organelles. Genetic studies in yeast have identified more than 40 vacuolar protein sorting (VPS) genes involved in vesicle transport to vacuoles. This gene is a member of the Sec1 domain family, and shows a high degree of sequence similarity to mouse, rat and yeast Vps45. The exact function of this gene is not known, but its high expression in peripheral



blood mononuclear cells suggests a role in trafficking proteins, including inflammatory mediators.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.