



# Recombinant Human Kelch repeat and BTB domain-containing protein 4 (KBTBD4)

<b>Product Code</b>	CSB-MP882119HU
<b>Abbreviation</b>	KBTBD4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NVX7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MESPEEPGAS MDENYFVNYT FKDRSHSGRV AQGIMKLCLE EELFADVTIS VEGREFQLHR LVLSAQSCFF RSMFTSNLKE AHNRVIVLQD VSESVFQLLV DYIYHGTVKL RAEELQEIYE VSDMYQLTSL FECSRFLAR TVQVGNCLQV MWLADRHSDP ELYTAAKHCA KTHLAQLQNT EEFLHLP HRL LIDIISDGVP CSQNPTEAIE AWINFNKEER EAFESLR TS LKEIGENVHI YLIGKESSRT HSLAVSLHCA EDDSI SVSGQ NSLCHQITAA CKHGGDLYVV GGSIPRRMWK CNNATVDWEW CAPLPRDLQ HTLVSVPGKD AIYSLGGKTL QDTLSNAVIY YRVGDNVWTE TTQLEVAVSG AAGANLNGII YLLGGEENDL DFFTKPSRLI QCFDTETDKC HVKPYVLPFA GRMHAAVHKD LVFIVAEGDS LVCYNPLLDS FTRLCLPEAW SSAPSLWKIA SCNGSIYVFR DRYKKGDANT YKLDPATSAV TVTRGIK VLL TNLQFVLA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	KBTBD4
<b>Protein Names</b>	Recommended name: Kelch repeat and BTB domain-containing protein 4 Alternative name(s): BTB and kelch domain-containing protein 4
<b>Expression Region</b>	1-518
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.