



Recombinant Human Required for meiotic nuclear division protein 1 homolog (RMND1)

Product Code	CSB-YP882123HU
Abbreviation	RMND1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NWS8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	HHILSKAH QCRRIGHLML KPLKEFENTT CSTLTIRQSL DLFLPDKTAS GLNKSQILEM NQKKSDTSM L SPLNAARCQD EKAHLPTMKS FGTHRRVTHK PNLLGSKWFI KILKRHFSSV STETFVPKQD FPQVKRPLKA SRTRQPSRTN LPVLSVNEDL MHCTAFATAD EYHLGNLSQD LASHGYVEVT SLPRDAANIL VMGVENSAKE GDPGTIFFFR EGAAVFWNVK DKTMKHVMKV LEKHEIQPYE IALVHWENEE LNYIKIEGQS KLHRGEIKLN SELDLDDAIL EKFAFSNALC LSVKLAIWEA SLDKFIESIQ SIPEALKAGK KVKLSHEEVM QKIGELFALR HRINLSSDFL ITPDFYWDRE NLEGLYDKTC QFLSIGRRVK VMNEKLQHCM ELTDLMRNHL NEKRALRLEW MIVILITIEV MFELGRVFF
Source	Yeast
Target Names	RMND1
Protein Names	Recommended name: Required for meiotic nuclear division protein 1 homolog
Expression Region	13-449
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.