



Recombinant Human Protein C20orf11 (C20orf11)

Product Code	CSB-EP882124HU-B
Abbreviation	C20orf11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NWU2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSYAEKPDEI TKDEWMEKLN NLHVQRADMN RLIMNYLVTE GFKEAAEKFR MESGIEPSVD LETLDERIKI REMILKGQIQ EAIALINSLH PELLDTNRYL YFHLQQQHLLI ELIRQRETEA ALEFAQTQLA EQGEESRECL TEMERTLALL AFDSPEESPFGDLLHTMQRQ KVVSEVNVQAV LDYENRESTP KLAKLLKLLL WAQNELDQKK VKYPKMTDLS KGVIEEPK
Source	E.coli
Target Names	GID8
Protein Names	Recommended name: Protein C20orf11 Alternative name(s): Two hybrid-associated protein 1 with RanBPM Short name= Twa1
Expression Region	1-228
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.