



# Recombinant Human FAST kinase domain-containing protein 2 (FASTKD2)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-BP882149HU  |
| <b>Abbreviation</b>      | FASTKD2   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q9NYY8  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Homo sapiens (Human)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MLTTLKPFGR VSVESKMNNK AGSFFWNLRQ FSTLVSTSRM MRLCCLGLCK<br>PKIVHSNWN LNNFHNRMQS TDIIRYLFQD AFIFKSDVGF QTKGISTLTA<br>LRIERLLYAK RLFVDSKQSL VPKVDSDEL KKVNLNHEVS NEDVLTKETK<br>PNRISRKLS EECNSLSDVL DAFSKAPTFP SSNYFTAMWT IAKRLSDDQK<br>RFEKRLMFSH PAFNQLCEHM MREAKIMQYK YLLFSLHAIV KLGIPQNTIL<br>VQTLRVTQE RINECDEICL SVLSTVLEAM EPCKNVHVL R TGFRILVDQQ<br>VWKIEDVFTL QVVMKCIGKD APIALKRKLE MKALRELD R SVLNSQHMFE<br>VLAAMNHRSL ILLDECSKVV LDNIHGCLP R IMINILQSCK DLQYHNLDLF<br>KGLADYVAAT FDIWKFRKVL FILILFENLG FRPVGLMDLF MKRIVEDPES<br>LNMKNLSIL HTYSSLNHVY KCQNKQFVE VMASALTGYL HTISSENLLD<br>AVYSFCLMNY FPLAPFNQLL QKDIISELLT SDDMKNAYKL HTLDTCLKLD<br>DTVYLRDIAL SLPQLPRELP SSHTNAKVAE VLSSLLGGEG HFSKDVHLP<br>NYHIDFEIRM DTNRNQVLPL SDVDTTSATD IQRVAVLCVS RSAYCLGSSH<br>PRGFLAMKMR HLNAMGFHVI LVNNWEMDKL EMEDAVTFLK TKIYSVEALP<br>VAAVNVQSTQ |
| <b>Source</b>            | Baculovirus   |
| <b>Target Names</b>      | FASTKD2   |
| <b>Protein Names</b>     | Recommended name: FAST kinase domain-containing protein 2   |
| <b>Expression Region</b> | 1-710   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Target Details</b>    | This gene encodes a protein that is localized in the mitochondrial inner compartment and that may play a role in mitochondrial apoptosis. Nonsense mutations have been reported to result in cytochrome c oxidase deficiency.   |

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.